

MODEL	T-	ONE	I-TECH	
VERSION	RE	НАВ	and anison	
 CONTENT Electrotherapy unit with g display and 4 channels. 4 connection cables (8 co electrodes). Self-adhesive square elect Self-adhesive rectangular Self-adhesive round elect Ionophoresis kit. Rechargeable Ni-MH battet Charger. User manual and electrod 		graphic colour onnectable ctrodes. r electrodes. trodes. ery.	Class IIa medical device	
manual. MAIN TECHNICAL FEATURES:				
 Four independent channels unit (up to 8 connectable electrodes); Constant current generator (can guarantee a constant current even when skin condition changes, the device automatically stops in the event of accidental electrodes disconnection); Compensated biphasic square wave: current quantity from positive pole to negative one is constant to avoid polarization effect; Monophasic wave generator (only for iontophoresis programs); 2+2 users programs; Backlight graphic colour display; User feedback: the device recognizes patient connection; Soft-touch easy to use keyboard; On - line guide on display for fast instructions; 4,8 V rechargeable battery pack; Personalized programs for man and woman; intensity adjustable from 1 to 120 mA; wave frequency from 0.2 to 1200 Hz; pulse width from 20 µs to 300 ms; time from 1 to 90 minutes; Medical device CE0068. 				
29 REHA PROGRAMS haematomas, oed			ar, Kotz, interferential, ionophoresis, microcurrent, na, Tens sequential, TENS Burst, TENS frequency/width prevention, atrophy)	
9 URO PROGRAMS		(stress incontinence, urgency incontinence, mixed incontinence)		

I.A.C.E.R. Srl

Via S. Pertini 24/A - 30030 Martellago (VE) - Italia / Italy

Tel.: (+39) 041/5401356 - Fax: (+39) 041/5402684 - Email: iacer@iacer.it - PEC: iacer@pec.it - Web: www.itechmedicaldivision.com Cod. Fisc. / P.IVA / Vat Number: IT00185480274 - Codice SDI: SUBM70N - R.E.A.: VE N. 120250 - M. VE001767 - Capitale Sociale / Share Capital: € 110.000,00 i.v.



16 TENS PROGRAMS	(TENS rapid, endorphinic, max values, inflammatory, backache/sciatic pain, sprains/bruises, vascularization,)		
20 NEMS PROGRAMS	(warming up, resistance, explosive strength, agonist/antagonist, muscle relaxant, deep massage, hypertrophy,)		
18 BEAUTY PROGRAMS	(firming up, toning up, definition, microlifting, lipolysis, modelling, capillarization,)		
22 2+2 USERS PROGRAMS	(TENS rapid, firming up, modelling, capillarization, resistance, basic strength, hypertrophy, microlifting, inflammatory,)		
12 FREE MEMORIES (MEM)	(free memories TENS, free memories NEMS)		
Dimensions: 150 x 90 x 35 mm			

SK329-01

I.A.C.E.R. Srl